

laboratory

glassware drying oven

ovens



 **labcon**
QUALITY ■ PRECISION ■ PERFORMANCE

labcon glassware drying oven

The glassware drying oven is made entirely from stainless steel 304, and is heated by sheathed incoloy electrical elements. The unit is fitted with 2 (two) large chromium plated wiremesh racks and features glass sliding doors front and back.

The temperature is adjustable by means of a Rheostat. A pilot light is fitted. The unit is supplied with a cable and plug.

technical data $\text{C}\epsilon$

MODEL	TEMPERATURE RANGE	CAPACITY LITRES	DIMENSIONS EXTERNAL			POWER RATING		WEIGHT Kg
			L	W	H	KW	VOLTAGE	
GDO 16	ambient to 120°C	160 l	1000	520	530	1000	220/240	42



Every effort has been taken to ensure the accuracy of this specification. However, due to continuous development and improvement the goods supplied may vary in detail from those described or illustrated.