

Precision  
Portable  
Incubator

laboratory



incubator

 labcon

QUALITY • PRECISION • PERFORMANCE

# labcon Precision Portable Incubator

LABCON Precision Portable Incubator Constant Temperature Storage unit have been designed and manufactured to the highest standards and include the following features:

- \* Temperature regulation to  $\pm 0.1^{\circ}\text{C}$
- \* Adjustable temperature setting from  $0^{\circ}$  to  $40.0^{\circ}\text{C}$ .
- \* Fan assisted air circulation
- \* Door switch with warning buzzer
- \* Three point temperature sensing
- \* 12V DC operation 75W maximum
- \* RS 232 interface for data capturing
- \* Temperature deviation alarm
- \* Operates over a  $+5^{\circ}\text{C}$  to  $+40^{\circ}\text{C}$  ambient temperature range
- \* Electronic cooling device. No gas refrigeration components.

#### OPTIONAL ITEMS:

- \* Digital display
- \* Cigarette lighter power cord
- \* Battery backup power supply
- \* Special shelves or racks as per customer specification.



- \* This unit is ideal for veterinarians that needs to travel with medicine or samples that needs to stay at a constant specified temperature.
- \* This Portable Incubator is also robust and assists in keeping samples at a specified temperature while other tests are being completed i.e. at a mobile testing station where the testing for the COVID-19 virus can be conducted without compromising the samples temperature or exposure to conditions that would compromise the samples.

## technical data

MODEL	PRODUCT CODE	DIMENSIONS						VOLUME (Lt)	SHELVES	POWER	WEIGHT (kg)
		EXTERNAL (mm)			INTERNAL (mm)						
		H	W	D	H	W	D				
PPI-01	ZPP01	610	400	430	444	275	325	40	4	220VAC/12VDC	6



**LABCON LABORATORY EQUIPMENT (PTY)LTD**  
 Reg. No: 2012/102814/07

P.O.Box 155, Maraisburg, 1700, South Africa  
 Phone: +27 11 762 7740 Fax: +27 11 762 8601  
 e-mail: [info@labcon.biz](mailto:info@labcon.biz) Website: [www.labcon.biz](http://www.labcon.biz)

Every effort has been taken to ensure the accuracy of this specification. However, due to continuous development and improvement the goods supplied may vary in detail from those described or illustrated.